



Dentsply Sirona

PMMA and Wax CAD/CAM Discs

MULTILAYER SHADED PMMA DISCS

Made from acrylic resins that are layered resulting in a highly esthetic provisional restoration. Multiple disc sizes available in all 16 Vita Shades and 1 bleach.

- Composed of 13 shade layers that blend seamlessly together
- Indicated for manufacturing temporary crowns and bridges
- Mills efficiently and polishes easily
- High flexural strength > 80 MPa
- Natural fluorescence and biocompatible

Size: 12mm

SHADE	PART #	SHADE	PART #
HW2 (BL2).....	8072612	B4	8071812
A1	8071012	C1	8071912
A2	8071112	C2	8072012
A3	8071212	C3	8072112
A3.5.....	8071312	C4	8072212
A4	8071412	D2	8072312
B1.....	8071512	D3	8072412
B2	8071612	D4	8072512
B3	8071712		

Size: 16mm

SHADE	PART #	SHADE	PART #
HW2 (BL2).....	8072616	B4	8071816
A1	8071016	C1	8071916
A2	8071116	C2	8072016
A3	8071216	C3	8072116
A3.5.....	8071316	C4	8072216
A4	8071416	D2	8072316
B1.....	8071516	D3	8072416
B2	8071616	D4	8072516
B3	8071716		

Size: 20mm

SHADE	PART #	SHADE	PART #
HW2 (BL2).....	8072620	B4	8071820
A1	8071020	C1	8071920
A2	8071120	C2	8072020
A3	8071220	C3	8072120
A3.5.....	8071320	C4	8072220
A4	8071420	D2	8072320
B1.....	8071520	D3	8072420
B2	8071620	D4	8072520
B3	8071720		

SHADED PMMA DISCS

These monochromatic PMMA discs deliver a strong, color stable provisional restoration. Multiple disc sizes available in all 16 Vita Shades, 1 bleach and 2 translucent shades and multiple disc sizes.

- High color stability
- Indicated for manufacturing temporary crowns and bridges
- Mills efficiently and polishes easily
- High flexural strength > 80 MPa
- Natural fluorescence and biocompatible

Size: 12mm

SHADE	PART #	SHADE	PART #
HW2 (BL2).....	8082612	B3	8081712
A1	8081012	B4	8081812
A1T	8083112	C1	8081912
A2	8081112	C2	8082012
A3	8081212	C3	8082112
A3.5.....	8081312	C4	8082212
A4	8081412	D2	8082312
B1.....	8081512	D3	8082412
B1T	8083212	D4	8082512
B2	8081612		

Dentsply Sirona
570 West College Avenue, York, PA 17401
1-800-243-1942

prosthetics.dentsply.com

©2016 Dentsply Sirona, Inc. All rights reserved.

DP-0000187 Rev. 1 (10/2016)



Size: 16mm

SHADE	PART #	SHADE	PART #
HW2 (BL2).....	8082616	B3	8081716
A1	8081016	B4	8081816
A1T	8083116	C1	8081916
A2	8081116	C2	8082016
A3	8081216	C3	8082116
A3.5.....	8081316	C4	8082216
A4	8081416	D2	8082316
B1.....	8081516	D3	8082416
B1T	8083216	D4	8082516
B2	8081616		

Size: 20mm

SHADE	PART #	SHADE	PART #
HW2 (BL2).....	8082620	B3	8081720
A1	8081020	B4	8081820
A1T	8083120	C1	8081920
A2	8081120	C2	8082020
A3	8081220	C3	8082120
A3.5.....	8081320	C4	8082220
A4	8081420	D2	8082320
B1.....	8081520	D3	8082420
B1T	8083220	D4	8082520
B2	8081620		

BURNOUT PMMA DISCS

Ideal for precise casting or pressing of metal and all-ceramic frameworks. This material burns out clean without leaving residue. Multiple disc sizes are available in red and clear.

- Material strength of >80 MPa results in greater precision in milling
- Excellent translucency
- Mills efficiently and polishes easily
- Burning residual < 0.05% (2h at 500°C)

SIZE	PART #
CLEAR	
12mm	8092712
16mm	8092716
20mm	8092720
25mm	8092725
RED	
12mm	8092812
16mm	8092816
20mm	8092820
25mm	8092825

WAX DISCS

This high quality Polyethylene wax is ideal for casting and pressing of metal and all-ceramic frameworks. It allows for easy carving when making adjustments and burns out completely without leaving any residue. Multiple disc sizes available in ivory and gray.

- Mills efficiently, resulting in a reduction of milling time and easy extraction
- Minimal CAM tool wear
- Burns out clean without residue
- Easy to carve and make adjustments

SIZE	PART #
IVORY	
16mm	8102916
20mm	8102920
25mm	8102925
GRAY	
16mm	8103016
20mm	8103020
25mm	8103025

